

Product Datasheet

CLK-IN-T3 (orb746306)

Description	CLK inhibitor T3 (T3) is a highly potent, selective, and cell-based stable CDC-like kinase (CLK) inh...
Note	For research use only
Formula	$C_{28}H_{30}N_6O_2$
MW	482.588
CAS Number	[2109805-56-1]
SMILES	<chem>O=C(C(C)(C)C1C=CC(C(NC2=CN3C=C(C4C=CN=CC=4)C=CC3=N2)=O)=CC=1)N1CCN(C)CC1</chem>
Expiration Date	12 months from date of receipt.