

## Product Datasheet

### TCPOBOP (orb746314)

<b>Description</b>	TCPOBOP is an agonist of the constitutive androstane receptor (CAR) (EC50 = 20 nM). TCPOBOP enhances...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{16}H_8Cl_4N_2O_2$
<b>MW</b>	402.052
<b>CAS Number</b>	[76150-91-9]
<b>SMILES</b>	<chem>C1C=CN=C(OC2=CC=C(OC3=NC=C(Cl)C=C3Cl)C=C2)C(Cl)=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.