



Product Datasheet

CTRP2 Antibody (orb74633)

Description	CTRP2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BXJ5
Expiration Date	12 months from date of receipt.



Western blot analysis of 3T3 (Balb) cell...



Immunohistochemical analysis of formalin...