

## Product Datasheet

### ML-179 (orb746345)

<b>Description</b>	ML-179 is a potent inverse agonist of liver receptor homolog-1 (LRH-1) (IC50 = 320 nM) with maximum ...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{25}F_3N_4O_2$
<b>MW</b>	422.44
<b>CAS Number</b>	[1883548-87-5]
<b>SMILES</b>	<chem>O=C(N1)N(C2CCCCC2)C(C=C1N(CC3)CCN3C4=CC=CC(C(F)(F)F)=C4)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.