

Product Datasheet

GSK-114 (orb746353)

Description	GSK-114 is a potent and selective selective TNNI3K inhibitor (TNNI3 IC ₅₀ = 25nM; B-RafV600E IC ₅₀ = 1000 nM). GSK-114 displays significant bias (40-fold) for TNNI3K over B-Raf, exceptional broad spectrum kinase selectivity, and adequate oral exposure to en
Note	For research use only
Formula	C ₁₉ H ₂₃ N ₅ O ₄ S
MW	417.48202252388
CAS Number	[1301761-96-5]
SMILES	S(C1C=CC(=C(C=1)NC1C2C=C(C(=CC=2N=CN=1)OC)OC)N(C)C)(NC)(=O)=O
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)