



Product Datasheet

CTRP7 Antibody (orb74642)

Description	CTRP7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8BVD7
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...



Immunohistochemical analysis of formalin...