

Product Datasheet

Alemtuzumab biosimilar antibody (orb746520)

Description	Alemtuzumab is a humanized monoclonal antibody specific to lymphocyte antigens. It is a recombinant ...
Conjugation	Unconjugated
Target	Alemtuzumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Purity	> 98%,Standard References
CAS Number	[216503-57-0]
Expiration Date	12 months from date of receipt.