

## Product Datasheet

### Panitumumab biosimilar antibody (orb746521)

<b>Description</b>	Panitumumab, formerly ABX-EGF, is a fully human monoclonal antibody specific to the epidermal growth...
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Panitumumab
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>CAS Number</b>	[339177-26-3]
<b>Expiration Date</b>	12 months from date of receipt.