

Product Datasheet

Efmoroctocog alfa biosimilar antibody (orb746522)

Description	Efmoroctocog alfa is a fully recombinant factor VIII-Fc fusion protein (rFVIII-Fc) with an extended h...
Conjugation	Unconjugated
Target	Efmoroctocog alfa
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Purity	> 98%,Standard References
CAS Number	[1270012-79-7]
Expiration Date	12 months from date of receipt.