

Product Datasheet

Natalizumab biosimilar antibody (orb746523)

Description	Natalizumab is a recombinant, humanized monoclonal antibody, binds to $\alpha4\beta1$ -integrin and blocks its...
Conjugation	Unconjugated
Target	Natalizumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Purity	> 98%, Standard References
CAS Number	[189261-10-7]
Expiration Date	12 months from date of receipt.