

Product Datasheet

Pertuzumab biosimilar antibody (orb746524)

Description	Pertuzumab, a humanized monoclonal antibody, is a HER2 dimerization inhibitor for the treatment of m...
Conjugation	Unconjugated
Target	Pertuzumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Purity	> 98%,Standard References
CAS Number	[380610-27-5]
Expiration Date	12 months from date of receipt.