

Product Datasheet

STING agonist compound 17 (orb746554)

Description	STING agonist compound 17 is a selective stimulator of interferon genes (STING) receptor agonist.
Note	For research use only
Formula	$C_{34}H_{40}N_{10}O_2$
Purity	> 98%, Standard References
MW	620.75
SMILES	<chem>C1C=CC=C2C=1N(/C(/N2CCCCN1C2=C(C=CC=C2)N(C)/C1=N\C(C1N(C)C)N=C(C)C=1)=O)=N\C(C1N(CC)N=C(C)C=1)=O)C</chem>
Expiration Date	12 months from date of receipt.