

## Product Datasheet

# Rabbit F(Ab')<sub>2</sub> Anti Goat IgG (H/L) (AP) ANTIBODY (orb746901)

<b>Description</b>	Rabbit polyclonal antibody against Rabbit F(Ab') <sub>2</sub> Anti Goat IgG (H/L) conjugated to AP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	AP
<b>Tested Applications</b>	DOT, ELISA, IHC-Fr, WB
<b>Target</b>	IgG
<b>Preservatives</b>	Purified IgG conjugated to Alkaline Phosphatase - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Alkaline Phosphatase - liquid
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.