

Product Datasheet

RABBIT F(ab')₂ ANTI GOAT IgG (H/L):Alk. Phos. (orb746901)

Description	Rabbit polyclonal antibody against Rabbit F(Ab') ₂ Anti Goat IgG (H/L) conjugated to AP
Species/Host	Rabbit
Reactivity	Goat
Conjugation	AP
Tested Applications	DOT, ELISA, IHC-Fr, WB
Target	IgG
Preservatives	Purified IgG conjugated to Alkaline Phosphatase - liquid
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)