

Product Datasheet

Rabbit F(Ab')₂ Anti Goat IgG (H/L) (FITC) ANTIBODY (orb746902)

Description	Rabbit polyclonal antibody against Rabbit F(Ab') ₂ Anti Goat IgG (H/L) conjugated to FITC
Species/Host	Rabbit
Reactivity	Goat
Conjugation	FITC
Tested Applications	ELISA, FC
Target	IgG
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.