

Product Datasheet

Rabbit F(Ab')₂ Anti Rat IgG (HRP) ANTIBODY (orb747149)

Description	Rabbit polyclonal antibody against Rabbit F(Ab') ₂ Anti Rat IgG conjugated to HRP
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	DOT, ELISA, IHC-Fr
Immunogen	Rat IgG.
Target	IgG
Preservatives	F(ab') ₂ fragment of IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Form/Appearance	F(ab') ₂ fragment of IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P20762 , P20759 , P20760 , P20761
Expiration Date	12 months from date of receipt.