



Product Datasheet

TWEAK Antibody (orb74715)

Description	TWEAK Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	043508
Expiration Date	12 months from date of receipt.



Western blot analysis of (A) 5 ng, (B) 2...



Immunohistochemical analysis of formalin...