



Product Datasheet

AIM Antibody (orb74720)

Description	AIM Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O43866
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji lysate at ...



Immunohistochemical analysis of formalin...