

Product Datasheet

HERPES SIMPLEX VIRUS 1 gG antibody (orb747242)

Description	Mouse monoclonal antibody to HERPES SIMPLEX VIRUS 1 gG
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, IH-Fr, IP, WB
Immunogen	HSV-1 strains HFEM and SC16.
Target	HERPES SIMPLEX VIRUS 1 gG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.