



Product Datasheet

BAG Antibody (orb74740)

Description	BAG Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q99933
Expiration Date	12 months from date of receipt.



Western blot analysis
of PC-3 cell lysat...



Immunohistochemical
analysis of formalin...