



Product Datasheet

Avian Influenza Hemagglutinin 4 Antibody (orb74757)

Description	Avian Influenza Hemagglutinin 4 Antibody
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q692M2
Expiration Date	12 months from date of receipt.



Western blot analysis of (A) 5 ng and (B...



The bar chart illustrates that 1ug/ml Av...