



Product Datasheet

sRANK Ligand Antibody (orb74770)

Description	sRANK Ligand Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O14788
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat liver tissu...



Immunohistochemical
analysis of formalin...