



Product Datasheet

VISA Antibody (orb74808)

Description	VISA Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q7Z434
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat brain tissu...



Immunohistochemical
analysis of formalin...