



Product Datasheet

Bub1 Antibody (orb74845)

| | |
|----------------------------|-------------------------------------------------------------------------------------------|
| Description | Bub1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, ICC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | 043683 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of A-20 cell lysat...



Immunohistochemical analysis of formalin...