

Product Datasheet

STREPTOLYSIN-O protein (orb748726)

Description	Recombinant STREPTOLYSIN-O protein
Reactivity	Bacteria
Tested Applications	ELISA
Target	STREPTOLYSIN-O
Preservatives	Recombinant purified Streptolysin O - liquid
Form/Appearance	Recombinant purified Streptolysin O - liquid
Concentration	Current, batch-specific concentration 3.04 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)