

Product Datasheet

Ferret IFN GAMMA protein (orb749115)

Description	Recombinant ferret IFN GAMMA protein
Reactivity	Ferret
Tested Applications	ELISA, WB
Target	IFN GAMMA
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	A4PIZ9
Expiration Date	6 months from date of receipt.