

Product Datasheet

Bovine TNF ALPHA protein (orb749117)

Description	Recombinant Bovine TNF ALPHA protein
Reactivity	Bovine
Tested Applications	ELISA, WB
Target	TNF ALPHA
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	A0A0U5EMB8
Expiration Date	6 months from date of receipt.