

Product Datasheet

Guinea pig TNF ALPHA protein (orb749118)

Description	Recombinant Guinea pig TNF ALPHA protein
Reactivity	Guinea pig
Tested Applications	ELISA, WB
Target	TNF ALPHA
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P51435
Expiration Date	6 months from date of receipt.