



Product Datasheet

FGF4 Antibody (orb74928)

Description	FGF4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08620
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH 3T3 cell ly...



Immunohistochemical analysis of formalin...