



Product Datasheet

Slitrk6 Antibody (orb74943)

Description	Slitrk6 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9H5Y7
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of SK-N-SH cell ly...