

Product Datasheet

p21 Antibody / CDKN1A (orb749443)

Description

This mAb recognizes a 21kDa protein, identified as the


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IF, IHC-P, WB



IHC: Formalin-fixed, paraffin-embedded h...

Immunogen


Recombinant full-length human protein was used as the immunogen for the p21 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the p21 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).




IHC: Formalin-fixed, paraffin-embedded h...

Note

For research use only

Application notes

Optimal dilution of the p21 antibody to be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min followed by cooling at RT for 20 min2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.



Immunofluorescent staining of human HeLa...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P38936

Hazard Information

This p21 antibody is available for research use only.

Dilution Range

Immunofluorescence: 1-2ug/ml,Western blot: 1-2ug/ml,Immunohistochemistry (FFPE): 2-4ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.