



Product Datasheet

CELA3B Antibody / Elastase 3B (orb749499)

Description	This mAb recognizes a protein of ~30kDa, identified as CELA3B (Chymotrypsin like elastase family mem...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A partial recombinant protein (aa 82-238) was used as the immunogen for the CELA3B antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the CELA3B antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the CELA3B antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P08861
Hazard Information	This CELA3B antibody is available for research use only.
Dilution Range	ELISA: 2-4ug/ml (order BSA/azide-free format), Western blot: 0.5-1ug/ml
Expiration Date	12 months from date of receipt.



Western blot of 1) partial recombinant p...



SDS-PAGE Analysis of Purified, BSA-Free ...