



## Product Datasheet

### CD15 Antibody / SSEA-1 (orb749512)

<b>Description</b>	CD15 plays a role in mediating phagocytosis, bactericidal activity, and chemotaxis. It is present on...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FACS, IF
<b>Immunogen</b>	Myelomonocytic leukemia cells were used as the immunogen for the FUT4 antibody.
<b>Preservatives</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
<b>Storage</b>	Store the FUT4 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
<b>Note</b>	For research use only
<b>Application notes</b>	Optimal dilution of the FUT4 antibody should be determined by the researcher.
<b>Formula</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
<b>Isotype</b>	Mouse IgM, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	PEG precipitation
<b>Uniprot ID</b>	<b>P22083</b>
<b>Hazard Information</b>	This FUT4 antibody is available for research use only.
<b>Dilution Range</b>	Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescent staining of human U937...



FACS analysis of human U937 cells using ...