


Product Datasheet

Interferon alpha 2 Antibody / IFNA2 (orb749587)

Description	Recognizes a protein of 16-27kDa, identified as human interferon-II) (IFN-II). Its epitope maps bet...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant human IFNA2 was used as the immunogen for the Interferon alpha 2 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the Interferon alpha 2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Interferon alpha 2 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P01563
Hazard Information	This Interferon alpha 2 antibody is available for research use only.
Dilution Range	ELISA (order the BSA-free format for coating)
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified,
BSA-Free ...