
Product Datasheet


Fas Antibody / CD95 (orb749603)

Description

mAb B-R18 specifically recognizes CD95, also known as

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, Functional Studies, IF, IHC-Fr
Immunogen	Recombinant human protein was used as the immunogen for the Fas antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the Fas antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Fas antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P25445
Hazard Information	This Fas antibody is available for research use only.
Dilution Range	Functional studies (order BSA/sodium azide-free format), Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (Frozen): 0.5-1ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.

NR 1



SDS-PAGE
Analysis of
Purified, BSA-
Free ...