



Product Datasheet

SATB1 Antibody (orb75002)

Description	SATB1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q01826
Expiration Date	12 months from date of receipt.



Western blot analysis of A20 cell lysate...



Immunohistochemical analysis of formalin...