



Product Datasheet

Napsin A antibody (Biotin) (orb750320)


Description	Mouse monoclonal antibody against Napsin A conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Amino acids 189-299 from the human protein were used as the immunogen for the Napsin A antibody.
Preservatives	0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Napsin A antibody at 2-8°C (up to one month) or aliquot and store at -20°C (longer term).
Note	For research use only
Application notes	Optimal dilution of the Napsin A antibody should be determined by the researcher.
Formula	0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	O96009
Hazard Information	This Napsin A antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 2-4ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human lung adenocar...



IHC staining of FFPE human lung adenocar...



Analysis of HuProt (TM) microarray conta...