
Product Datasheet

Human Coagulation Factor III / Tissue Factor Protein (orb750336)

Description

Human Tissue Factor, Fc Tag is

(HEK293). It contains AA Ser 33 - Glu ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Coagulation Factor III

Preservatives

50 mM Tris, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.2 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_001984.1

Purity

95%

Source

Human Tissue Factor, Fc Tag (orb750336) is expressed from human 293 cells (HEK293). It contains AA Ser 33 - Glu 251 (Accession # P13726-1).

MW

51.2 kDa

Uniprot ID[P13726](#)**NCBI**[NP_001984.1](#)**Expiration Date**

6 months from date of receipt.