



Product Datasheet

Human PTK7 / CCK4 Protein (orb750342)

Description	Human PTK7, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 31 - Thr 704 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PTK7
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 76.5 kDa. The protein migrates as 90-115 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002812.2
Purity	95%
Source	Human PTK7, His Tag (orb750342) is expressed from human 293 cells (HEK293). It contains AA Ala 31 - Thr 704 (Accession # Q13308-1).
MW	76.5 kDa
Uniprot ID	Q13308
NCBI	NP_002812.2
Expiration Date	6 months from date of receipt.



Flow
cytometric
analysis of
Human PTK7
p...



SDS-PAGE
analysis of
Human PTK7
protein.