

Product Datasheet

Human CD45 / PTPRC protein (orb750349)

Description	Human CD45, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Lys 577 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD45
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 62.8 kDa. The protein migrates as 120-170 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002829.3
Purity	90%
Source	Human CD45, His Tag (orb750349) is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Lys 577 (Accession # P08575-3).
MW	62.8 kDa
Uniprot ID	P08575
NCBI	NP_002829.3
Expiration Date	6 months from date of receipt.