

Product Datasheet

Human ZNRF3 protein (orb750350)

Description Human ZNRF3, His Tag is expressed from human 293 cells (HEK293). It contains AA Lys 56 - Met 219 (Ac...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target ZNRF3

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 19.7 kDa. The protein migrates as 22-25 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001193927.1

Purity 95%

Source Human ZNRF3, His Tag (orb750350) is expressed from human 293 cells (HEK293). It contains AA Lys 56 - Met 219 (Accession # Q9ULT6-1).

MW 19.7 kDa

Uniprot ID [Q9ULT6](#)

NCBI [NP_001193927.1](#)

Expiration Date 6 months from date of receipt.