

Product Datasheet

Human LGR4 / GPR48 protein (orb750351)

Description Human LGR4, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 25 -Thr 544 (Acce...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target LGR4

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.8 kDa. The protein migrates as 70-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_060960.2

Purity 85%

Source Human LGR4, His Tag (orb750351) is expressed from human 293 cells (HEK293). It contains AA Ala 25 -Thr 544 (Accession # Q9BXB1-1).

MW 58.8 kDa

Uniprot ID [Q9BXB1](#)

NCBI [NP_060960.2](#)

Expiration Date 6 months from date of receipt.