

Product Datasheet

Human LGR5 / GPR49 protein (orb750352)

Description

Human LGR5, His Tag is expressed from

Gly 22 - Arg 561 (Acc...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LGR5

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 62.0 kDa. The protein migrates as 65-75 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_003658.1

Purity

90%

Source

Human LGR5, His Tag (orb750352) is expressed from human 293 cells (HEK293). It contains AA Gly 22 - Arg 561 (Accession # O75473-1).

MW

62.0 kDa

Uniprot ID
O75473
NCBI
NP_003658.1
Expiration Date

6 months from date of receipt.