

Product Datasheet

Human PLGF / PGF (19-149) Protein (orb750367)

Description	Human PIGF-1 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 19 - Arg...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PLGF
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 41.2 kDa. The protein migrates as 50-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001193941.1
Purity	95%
Source	Human PLGF (19-149), Fc Tag (orb750367) is expressed from human 293 cells (HEK293). It contains AA Leu 19 - Arg 149 (Accession # P49763-2).
MW	41.2 kDa
Uniprot ID	P49763
NCBI	NP_001193941.1
Expiration Date	6 months from date of receipt.