

Product Datasheet

Human NOTCH1 Protein (orb750370)

Description	Human NOTCH1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Gln 526 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	NOTCH1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 80.0 kDa. The protein migrates as 90-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_060087.3
Purity	95%
Source	Human NOTCH1, Fc Tag (orb750370) is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Gln 526 (Accession # P46531-1).
MW	80.0 kDa
Uniprot ID	P46531
NCBI	NP_060087.3
Expiration Date	6 months from date of receipt.