



Product Datasheet

TRIM30 Antibody (orb75038)

Description	TRIM30 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P15533
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse spleen ti...



Immunohistochemical analysis of formalin...