



Product Datasheet

CDCA8 Antibody (orb75041)

Description	CDCA8 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q53HL2
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat kidney lysa...



Confocal immunofluorescence analysis of ...