



Product Datasheet

POFUT1 Antibody (orb75047)

Description	POFUT1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9H488
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...



Immunohistochemical analysis of formalin...