

Product Datasheet

Morc3 antibody (orb750550)

Description

Morc3 antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Morc3 whole rabbit serum was prepared by repeated immunizations with a mouse Morc3 recombinant protein.

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Liquid (sterile filtered)

Concentration

80 mg/mL

Storage

Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Morc3 Antibody is tested for use in Western Blot and suitable for ChIP and IF. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 110.8 kDa in size corresponding Morc3 by western blotting in the appropriate cell lysate or extract.

Isotype

Antiserum

Clonality

Polyclonal

Purity


Morc3 antibody was prepared from monospecific antiserum by delipidation and defibrination. The antibody is specific for mouse Morc3 in expressed cell lysates. Cross reactivity to other Morc proteins has not been determined.

Uniprot ID
[Q6A0C2](#)
NCBI
[NP_001038994.2](#)
Dilution Range

ELISA: 1:10,000-1:50,000, WB: 1:1000-1:5000

Expiration Date

12 months from date of receipt.



Western Blot of Rabbit Anti-Morc3 antibo...