
Product Datasheet

Rabbit IgG F(ab')₂ (orb750722)

DescriptionRabbit IgG F(ab')₂**Species/Host**

Goat

Reactivity

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, WB

ImmunogenRabbit IgG F(ab')₂ fragment**Preservatives**

None

Form/Appearance

Lyophilized

Concentration

31.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Rabbit IgG F(ab')₂ antibody is suitable for use in ELISA, immunohistochemistry, and western blot. Specific conditions for reactivity should be optimized by the end user.

Isotype

Antiserum

Clonality

Polyclonal

Purity

This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-goat serum, Rabbit IgG, Rabbit IgG F(ab')₂ and Rabbit Serum. No reaction was observed against Rabbit IgG F(c).

Dilution Range

ELISA: 1:20,000 - 1:100,000, IHC: 1:1,000 - 1:5,000, WB: 1:2,000 - 1:10,000

Expiration Date

12 months from date of receipt.